PENIS EXTENDER

EAN-Code: 8714273926005

Markenname: Pumped by Shots

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Einführung

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Instructions for use

Introduction

The Penis Extender is designed to help you achieve your desired results safely and comfortably. Please read the following instructions carefully before use to ensure proper application and optimal performance.

Components

- Penis Extender device
- Adjustable straps
- Comfort pads

Instructions for Use

1. Preparation

- Ensure that the Penis Extender and all its components are clean and dry before use.
- Find a comfortable and private space to apply the device.

2. Assembly

- Attach the adjustable straps to the main body of the Penis Extender.
- Place the comfort pads on the areas that will come into contact with your skin to prevent discomfort.

3. Application

- Gently slide your penis into the Penis Extender, ensuring that it is flaccid during application.
- Adjust the base ring so that it sits comfortably at the base of your penis.
- Secure the adjustable straps around your shaft, ensuring they are snug but not too tight to avoid restricting blood flow.

4. Adjustment

- Gradually increase the tension using the adjustment mechanism. Start with a light tension and slowly increase over time as your comfort level allows.
- It is recommended to use the Penis Extender for short periods initially, gradually increasing the duration as you become accustomed to the device.

5. Duration and Frequency

- Use the Penis Extender for a few hours a day, as per your comfort level.
- Consistency is key, so aim to use the device regularly for optimal results.

6. After Use

- Carefully remove the Penis Extender by loosening the straps and sliding it off.
- Clean the device with mild soap and water, then dry thoroughly before storing.

7. Safety and Maintenance

- Regularly inspect the device for any signs of wear or damage and replace parts if necessary.
- Discontinue use immediately if you experience any pain or discomfort, and consult a healthcare professional if necessary.

Important Notes

- The Penis Extender is intended for personal use only and should not be shared with others.
- Keep the device out of reach of children.
- For external use only. Do not use on broken or irritated skin.

By following these instructions, you can safely and effectively use your Penis Extender to work towards your enhancement goals.

Washing instructions

To ensure the longevity and hygiene of your Penis Extender, please follow these cleaning instructions:

1. Before First Use: Wash the product with warm water and mild soap to remove any manufacturing residues. Rinse thoroughly and dry with a clean towel.
2. After Each Use:
 - Carefully detach the Penis Extender from any accessories or attachments.
 - Rinse the product under warm running water to remove any bodily fluids, lubricants, or other residues.
 - Use a gentle, unscented soap to clean the surface of the extender. Avoid using harsh chemicals or abrasive materials that may damage the product.
 - Rinse thoroughly to ensure all soap is removed.
3. Drying: After washing, pat the Penis Extender dry with a clean towel. Alternatively, allow it to air dry completely before storing.
4. Storage: Store your Penis Extender in a cool, dry place, away from direct sunlight and extreme temperatures.
5. Additional Care: For deeper cleaning, you may soak the product in a solution of warm water and mild soap for a few minutes, then rinse thoroughly.

By following these instructions, you will help maintain the quality and hygiene of your Penis Extender for continued use.

Ingredients

-

Warnings

Before using the Penis Extender, please read and understand the following warnings and precautions:

1. Consult a Healthcare Professional: It is strongly recommended that you consult with a qualified healthcare professional before using this product, especially if you have any pre-existing medical conditions, are taking medications, or have concerns about your sexual health.
2. Potential Discomfort: Some users may experience discomfort, pain, or irritation during or after use. If you experience severe pain or discomfort, discontinue use immediately and consult a healthcare professional.
3. Use with Caution: Improper use of the Penis Extender can lead to injury, including but not limited to bruising, swelling, or tissue damage. Always follow the manufacturer's instructions for safe and effective use.
4. Not for Everyone: This product is not recommended for individuals with certain medical conditions, such as Peyronie's disease, erectile dysfunction, or any other genital or penile disorders.
5. Monitor Your Progress: Regularly monitor your body's response to the product. If you notice any abnormal changes or reactions, cease use immediately and seek medical advice.

6. Age Restrictions: This product is intended for adult use only. Do not use if you are under the age of 18.

7. Storage and Maintenance: Ensure proper cleaning and storage of the product to prevent hygiene-related issues.

By using the Penis Extender, you acknowledge that you have read, understood, and accepted these warnings. Your safety and well-being are our top priority.

Disposal instructions

When it comes time to dispose of your Penis Extender, please follow these guidelines to ensure proper disposal and recycling:

1. **Disassemble the Product:** If your Penis Extender consists of multiple parts, carefully disassemble them to facilitate proper disposal.
2. **Check Local Regulations:** Before disposal, check your local waste management regulations to determine if there are specific instructions for disposing of personal items.
3. **Non-Recyclable Components:** Any non-recyclable components, such as plastic or silicone parts, should be placed in your regular household waste. Ensure they are securely wrapped to prevent any harm to waste management personnel.
4. **Recyclable Materials:** If your Penis Extender includes any recyclable materials (such as metal components), please separate these and place them in the appropriate recycling bin according to your local recycling guidelines.
5. **Hygiene Considerations:** For hygiene purposes, it is recommended to clean the product thoroughly before disposal. This helps prevent any potential contamination.
6. **Environmental Responsibility:** If you are unsure about how to dispose of the product responsibly, consider contacting your local waste disposal facility for guidance on environmentally friendly disposal options.

By following these disposal instructions, you help contribute to a cleaner environment while ensuring safe and responsible waste management. Thank you for your cooperation!

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declare under our sole responsibility for the product(s):

Equipment:	Penis Extender
Brand name:	Pumped by Shots
Model No.:	PMP029BLK
EAN code:	8714273926005

that all components of the above-mentioned product(s) are compliant with the following legal requirements:

REACH

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the registration, evaluation, authorisation and restriction of chemicals pursuant to article 33 (1) (**SVHC**) and **Annex XVII** of (EC) No 1907/2006

GPSR

General Product Safety Regulation (EU) 2023/988 (GPSR)

Cal Prop 65

California Proposition 65 (Safe Drinking Water and Toxic Enforcement Act of 1986), with reference to lead, phthalates and cadmium

PFAS Statement

We hereby declare that we do not intentionally add perfluoroalkyl and polyfluoroalkyl substances (PFAS) or any associated chemicals as a raw material or processing aid in the manufacturing of our products and their packaging.

Signed:

A handwritten signature in black ink, appearing to read "Hein Schouten".

Name: Hein Schouten

Position: CEO

Date: 16th of October, 2023